

# Use of Fatty Acids for the Identification of Phytopathogenic Bacteria

MARGARET A. ROY, Microbial ID, Inc., Newark, DE

Whole-cell fatty acid (FA) profiles have been used in bacterial classification for over 20 years and have become increasingly important in bacterial identification. Any microbiologist with access to a gas chromatograph (GC) can perform FA analysis.

Under standardized conditions, bacterial FAs are stable and reproducible within taxa (though they are modified in predictable fashion with environmental changes.) The FA profiles are highly diagnostic for a given species regardless of geographical source, though for some groups of pathogens there are subgroups related to host range. Fatty acid-based groupings have so far agreed with DNA homology, where DNA homology data are available.

The FAs between nine and 20 carbon lengths are most useful in identification of plant-associated bacteria. Presently, over 200 FAs are used for identification purposes, and these include saturated, unsaturated, hydroxy, cyclopropane, iso, and anteiso FAs (3).

**Methods.** For pure cultures, approximately 40 mg of wet cells is placed in a screw-capped test tube, and the FAs in whole cells are saponified, methylated, extracted into hexane/ether, and chromatographed (2,3).

It is possible to extract bacterial FAs directly from infected plant tissue for identification to the genus level and often to the species level, taking care to exclude as much plant lipid as possible. Since the bacteria grew in plant tissue, their FA profiles will be slightly different from their counterparts grown under controlled laboratory conditions. Improved methods for this extraction technique are being developed.

Any GC can be used for identification of bacteria, but capillary column GC (5% phenyl methyl silicone) gives the best resolution of peaks. Without a computer, chromatograms of unknowns can be compared with those of type or known strains by placing one chromatogram over the other and holding both up to a light source to see if the patterns are similar or different. This will distinguish genera—for example, *Xanthomonas* from *Pseudomonas*—but may not permit resolution at the species or pathovar level. Additionally, the FAs responsible for each peak are not identified or quantified. The individual component FAs can be identified by comparisons with FA standards purchased from fine chemical supply houses. FA names, peak relative retention time, and overall profiles for known strains published in the literature can also be used as references. For most genera, qualitative differences in FAs can be used for identification. For best results at the species or subspecific levels, the ratios between FAs that are most stable (i.e., the straight chain, iso, and anteiso FAs) should be calculated.

Over 17,000 samples per year can be analyzed on an automated GC system with an automatic sampler, an integrator, and a computer (Microbial ID, Inc., Newark, DE). A computer program calculates the peak areas and compares the unknown FA profiles with those of known reference strains in the data base. A similarity index expresses how nearly the profile of the unknown bacterium fits the closest species match in the data base. At present, the data bases contain FA profiles of over 300 species of aerobic bacteria and over 100 species of anaerobic bacteria. The software allows the user to expand or create data bases by adding FA profiles of additional species.

**Examples of fatty acids in phytopathogenic bacteria.** Within *Agrobacterium tumefaciens*, the three biotypes can be distinguished primarily by relative amounts of cyclopropane FAs. The tumor-inducing plasmid does not affect whole-cell FA profiles.

*Clavibacter michiganense* subsp. *michiganense* has a simple profile characterized by iso and anteiso FAs. To delineate subspecies, ratios between major peak areas should be calculated.

Most members of the Enterobacteriaceae contain the FA 14:0 3OH. *Erwinia carotovora* subsp. *atroseptica* is distinguished from *E. c.* subsp. *carotovora* by the ratios of 12:0 to 14:0, of 16:0 to 12:0, and of 16:1 to 18:1 (1).

The fluorescent pseudomonads seldom contain more than a trace amount of 14:0 3OH (a major peak in *Pseudomonas cepacia* and related nonfluorescent bacteria) but do contain 10:0 3OH. Within *P. syringae*, it is possible to distinguish pathovars such as, for example, *P. s.* pv. *tomato* and *P. s.* pv. *papulans* from other pathovars of *P. syringae* on the basis of presence (others) or absence (*P. s.* pv. *tomato* and *P. s.* pv. *papulans*) of 17:0 cyclopropane.

In the genus *Xanthomonas*, 11:0 ISO, 11:0 ISO 3OH, 13:0 ISO 3OH, 15:0 ISO, and 15:0 ANTEISO are important components of a very complex profile, while calculation of ratios between 15:0 ISO and 15:0 ANTEISO is necessary to determine the species. *X. campestris* pathovars are by definition indistinguishable using biochemical tests, yet xanthomonads can often be discriminated at the pathovar level using ratios of FAs.

#### LITERATURE CITED

1. De Boer, S. H., and Sasser, M. 1986. Differentiation of *Erwinia carotovora* ssp. *carotovora* and *E. carotovora* ssp. *atroseptica* on the basis of cellular fatty acid composition. *Can. J. Microbiol.* 32:796-800.
2. Miller, L., and Berger, T. 1985. Bacteria identification by gas chromatography of whole cell fatty acids. Hewlett-Packard Application Note 228-41. 8 pp.
3. Sasser, M. 1988. Identification of bacteria through fatty acid analysis. In: *Methods in Phytobacteriology*. Z. Klement, K. Rudolf, and D. C. Sands, eds. Akademiai Kiado, Budapest. In press.

## Salute to APS Sustaining Associates

This section is designed to help APS members understand more about APS Sustaining Associates. Information was supplied by company representatives. Each month different companies will be featured. A complete listing appears in each issue of *Phytopathology*.

**W. R. Landis Associates, Inc., Contact: W. R. Landis, President, P.O. Box 5126, Valdosta, GA 31603-5126; 912/247-6472.** W. R. Landis is a research and development management company that specializes in agricultural products. We have the capability of guiding products through their entire research and development process, with emphasis on the requirements of the U.S. Environmental Protection Agency (EPA) and state regulatory agencies. Because we serve clients from the smallest formulators to the largest multinational giants, our regulatory expertise with the EPA and individual states is broad-based. We keep abreast of and can conduct all federal and state requirements under Good Laboratory Practices (GLP) guidelines. These requirements include formulation development and evaluation, toxicology, environmental chemistry, and nontarget organism studies. Landis takes a project management approach to research and development and has the disciplines of plant pathology, entomology, agronomy, nematology, and chemistry on staff to bring technical expertise to each project. In addition, the company has expanded its in-house capabilities by acquiring Southern Agricultural Research, Inc., with two field research laboratories that offer full-spectrum contract field research in Georgia, Alabama, and Florida, and because of the demands of the midwest and west, we have recently acquired Northwest Agricultural Research, Inc., with the capability of conducting research in Idaho, Montana, Oregon, and Washington. We are also in the final stages of establishing facilities in Ohio and California.

**Eli Lilly & Company, Greenfield Lab, Contact: Bill Wright, P.O. Box 708, Greenfield, IN 46140; 317/467-4118.**

**Merck Company Inc., Contact: Richard A. Pence, P.O. Box 2000, Rahway, NJ 07065; 201/750-8605.**

**Mobay Corporation, Contact: K. A. Noegel, Research Fungicide Product Manager, Agricultural Chemicals Division, Kansas City, MO 64120; 816/242-2752.** Mobay is a basic manufacturer of a diversified line of agricultural chemicals and specialty chemical products. The Agricultural Chemicals Division is a leader in the development and marketing of ergosterol biosynthesis inhibiting fungicides in the United States. One of these products is Bayleton Systemic Fungicide, which is registered for the control of specific diseases in agricultural and specialty crops. Mobay also manufactures and distributes the fungicide Dyrene, which is used for the control of certain diseases in vegetable crops and turfgrasses.

**Monsanto Agricultural Company, Contact: Brad Castanho, Manager, 700 Chesterfield Villlage Parkway (GG3A), Chesterfield, MO 63198; 314/537-6663.** Monsanto Agricultural Company is an operating unit of Monsanto Company, headquartered in St. Louis, MO. It is a recognized world leader in the research, development, manufacturing, and marketing of crop protection chemicals. The company also produces a variety of animal feed supplements and preservatives and owns two seed companies that produce high-quality soybean and wheat seed. Total 1986 sales were more than \$1.1 billion, with annual research and development expenditures exceeding \$100 million. Areas of new and product research include chemical weed disease and insect control, and an increasing emphasis on applying the tools of biotechnology to crop improvement and protection and to animal nutrition and health.

---

## 1988 Advertisers Index

Academic Press .....	
LI-COR .....	373
Springer-Verlag New York, Inc. ....	
U.S. Department of Agriculture .....	